



Culique

Thermal Frequency Hair
Growth Treatment

Culique Hat

QUICK START GUIDE



Culique
Infrared
Hat



Culique
Controller



Power &
USB -C
Cable

How to Use

- 1 Unbox the Components**
Open the package and carefully take out all items.
Make sure all components are included and ready for setup.



- 2 Select the Right Plug**
Choose the appropriate plug adapter for your region.
Attach it securely to the power supply unit.



- 3 Connect Controller to Power**
Plug the power adapter into a power source.
Ensure the controller is receiving power.



- 4 Connect USB-C to Controller**
Insert the USB-C cable into the controller port.
Make sure the connection is firm and secure.



- 4 Connect Cable to Hat**
Attach the other end of the cable to the Culique Hat.
Check that the connection is properly aligned.



- 5 Culique Hat – Ready to Use**
Once everything is connected, the device is ready.
Turn on the controller and begin your session.

